# EZ 5500

---

## Bedienungs- anleitung/Garantie

Gebbruksaanwijzing  
Mode d'emploi  
Instrucciones de servicio  
Manual de instruções  
Istruzioni per l'uso  
Bruksanvisning  
Instruction Manual  
Instrukcja obsługi/Gwarancja  
Návod k použití  
Használati utasítás  
Руководство по эксплуатации

## Elektrische Zahnbürste

Elektrische tandenborstel  
Brosse à dents électrique  
Cepillo de dientes eléctrico  
Escova de lavar dentes eléctrica  
Spazzolino da denti elettrico  
Elektrisk tannbørste  
Electric Toothbrush  
Elektryczna szczoteczka do zębów  
Elektrický zubní kartáček  
Elektromos fogkefe  
Электрическая зубная щётка

# AEG

PERFECT IN FORM AND FUNCTION

## 2 Deutsch

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---

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## Safety Information

---

- The device is used for dental and oral care in adults.
- Children from the age of seven may use the device if supervised by an adult.

**WARNING:**

This machine is not suitable for children under the age of 3 as it contains small parts which can be broken off and swallowed.

## Installing the Battery

---

The battery compartment on the bottom of the device can be unlocked with a coin or similar item. Then you can remove the lower section of the housing. Insert the batteries (2x 1.5 V AA). Please ensure the correct polarity.

### Batteries (not supplied)

If the device is not used for a long period of time, remove the batteries to avoid the spilling of battery acid.

Different battery types or new and used batteries must not be used together.

**CAUTION:**

Batteries are not to be disposed of together with domestic wastes. Dispose of the used batteries at the relevant collection bin or at your dealer's.

**WARNING:**

Do not expose the batteries to high temperatures or direct sunlight. Never throw batteries into the fire. There is a **danger of explosion!**

## Instructions for Use

---

Insert a brush diagonally on to the drive rod and turn it until it locks in place. The bristles must be pointing in the direction of the burls. The brushes are fitted with rings of different colours for allocation to different people.

**i** NOTE:

The toothbrush does not start to operate immediately when you press the ON/OFF switch, but has to be pressed against the teeth at the same time.

Put some toothpaste on the brushes and switch the device on. Move the head of the brush carefully from tooth to tooth.

The brush should follow both the edge of the gums and the spaces between the teeth. It is not necessary to press firmly.

## Cleaning

---

In order to release the brush, turn the head of the brush in a clockwise direction. Now you can remove it completely. Please do not use any cleaning agents.

Net weight:.....0,12 kg

## Environmental Protection

---

When you wish to dispose of the device, please do not put it in the standard domestic waste. Instead you should take it to an official collection point for recycling. In this way you help to protect the environment.

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