**Untere Abb.:** Batteriewechsel  
**Fig. below:** Battery replacement  
**Afb. hieronder:** Vervangen van de batterij  
**Fig. ci-dessous :** Remplacement de la pile  
**Fig. a continuación:** Sustitución de la batería  
**Fig. sotto:** Sostituzione della batteria  
**Ilustracja poniżej:** Wymiana baterii  
**Alábbi abra:** Az elem cseréje  
**Мал. niche:** Замена батареек  
**Рис.:** Замена батареек

## 1 Deutsch

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### User manual

#### Please read this manual before use!

Thank you for your choice of our clinical digital thermometer FT 3617 (FDTH-VO-3). Famidoc Technology specializes in the research, development and selling of health care devices and baby care products. Our digital thermometers have been sold to the market of Europe and America for years with good reputation. After the research and investigation of the market response, it is found that 90% of the inaccuracy and function failure are resulted from improper usage. Thus, please read this manual carefully before using this product. Please keep your manual to avoid trouble for your usage.

### Warning, notice and suggestion

Clean and sanitize the probe of thermometer by 75% medical alcohol before and after the measurement each time.	Without prior disinfection, the use of thermometers may communicate pathogenic bacteria.
Do not allow children to measure body temperature by themselves. Measurement must be accompanied by an adult.	The children own operation may hurt themselves or cause inaccurate measurement.
This thermometer is only intended for oral, underarm and rectal temperature measurement. Don't use it for other intended use.	The unintended usage may hurt or mislead you or the patient.
The measurement can only be used for health monitoring and reference. Please don't conduct diagnose and treatment based on the measurement by yourself.	The thermometer can't replace the doctor's diagnosis. If you find the body temperature is abnormal, please contact your doctor.
The thermometer is not a toy, store it in a safe place out of children's reach. Don't disassemble, warp, grip or drop the thermometer.	Rude usage can damage the thermometer; Playing improperly by children may it or even swallow the parts such as the battery and steel cap. Please send them to hospital as soon as possible if that happens.
Don't throw the thermometer or the battery into fire.	The battery is likely to explode in fire!
Don't insert the metal probe of the thermometer into the plug of the power supply.	This false operation may lead to electric shock.
Don't use the sharp tool, scissors etc. to scissor, cut, pierce, scratch the measuring probe.	When you use again, the measuring probe parts may be pealed off and swallowed by the bite and pull.

### Notice

**!** Don't take temperature within 30 minutes after sports, bath and swimming.

## 6 English

In order to ensure the accuracy of measurement, please don't measure body temperature in strong electromagnetic interference environment (such as short wave, high frequency equipment operation place).

The improper use and maintenance to thermometer can be harmful to the human body, or leads to partial loss of function because the acceleration component aging.

Do not expose the thermometer directly to the sun, too wet, presence of flammable gases, dust and extreme variable temperature environment.

Don't take temperature after drinking hot water.

Don't take oral temperature when you have dental ulcer.

Don't take rectal temperature when you have hemorrhoids.

Do not measuring underarm if there are large area of scar. Scar connective tissue will affect the real temperature of the reaction, resulting in low temperature.

Don't take temperature immediately when you take medicine which may increase the body temperature.

The thermometer must be stored according to the technical specifications. Do not measure while talking on the phone.

### Recommendations

Don't try other method to disinfect the thermometer besides the method this manual prescribed.  
 Stop using the thermometer when there is any allergy or discomfort.  
 Take probe cover when taking oral temperature for hygienic reason.  
 Normally, female temperature during period is higher than usual by 0.3-0.5 °C.

### Intended Use

This clinical digital thermometer is determined for measuring the human body temperature in the mouth (oral), the armpit (axillary) or the anus (rectal). This thermometer is suitable for people at different ages, including infants, children and adults. In order to ensure the accuracy and safety of the measurement, adults should take temperature for children.

### Side effect

Nausea and discomfort may appear among the sensitive group when taking oral and rectal temperature. Both are normal phenomenon and can be relieved after rest.

### Feature

- Delicate ergonomic design, easy to use.
  - Last memory, convenient to check last measurement.
  - 1 minute readout the temperature accurately.
  - Convenient and accurate
  - Microcomputer control with buzz indication function, easier use.
  - Automatic switch off
- If forgetting to turn off the thermometer, it would enters power saving condition to prolong the service life of battery.
6. Environment—friendly
- Non glassy nor mercury, avoid pollution to the environment and skin. The probe and body materials conform to requirement of biocompatibility to ensure safety.

### Product Description

**Structure:**  
 1. Measurement tip / 2. Unit Body / 3. LCD / 4. Power button / 5. Battery cap

#### Battery install and change

Hold the body tightly with one hand and hold the battery cap with the other hand. Remove the battery cap as the figure "Battery replacement" shows. Then use a small tool to take out the spent AG3 battery and replace a new one. Don't install the battery in a reverse direction. Put down the side with "+" signal.  
 Make sure the AG3 battery is installed properly, then put on the battery cap as before.  
**NOTICE: Spent batteries must be disposed of according to the native law.**

### Operation instruction

**Display information**

6. Temperature / 7. Unit of Temperature (°C) / 8. Low battery symbol

#### Instrutions

- Press the [Power] button, then the LCD will show (18.88)° for 2 seconds, every time press the button, you will hear a "bi".
- Then the LCD will display last measurement for 2 seconds (such as 36.6° or 100.0°). If the power is too low and last measurement is effective, then the LCD will flashes to indicate low power (such as 36.6°). If the power is too low while last measurement is not effective, then LCD wont show the low power indication signal (such as 100.0°).
- When the memory display ends, LCD will show the automatically test temperature about 37.0° one second. At the same time, the temperature unit °C will flash.
- When LCD finishes self temperature test display, the thermometer is ready to work. When the measurement keep still for 16 seconds, the °C signal will stop flashing followed with 10 times Bi-Bi alarm sound to indicate the finish of measurement. Don't take out the thermometer within 3-5 minutes in order to get accurate measurement result.
- Automatic switch off  
 If forgotten to be switched off after the measurement, it will automatically shut off within 10 minutes.

#### Measurement method

- In the mouth (oral)

## 7 English

Position the thermometer in one of the two pockets under the tongue, to the left or right at the root of the tongue. The measuring probe must be in good contact with the tissue. Close your mouth and breath evenly through the nose to prevent the measurement from being influenced by inhaled/ exhaled air.

To consider the healthy requirements, please take the probe cover when taking the temperature. We don't sell probe cover, please buy them from professional thermometer market with CE certificate and comply with biocompatibility requirements.

- In the anus (rectal)  
 This is the most reliable measurement method, and is especially suitable for infants and children. Carefully insert the measuring probe of the thermometer

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Famidoc Technology Co., Ltd.
Shanghai International Holding Corp GmbH (Europe)