MS 5504

Bedienungsanleitung/Garantie

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General Safety Instructions

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The appliance is designed exclusively for private use and for the envisaged purpose. This appliance is not fit for commercial use. Do not use it outdoors. Keep it away from sources of heat, direct sunlight, humidity (never dip it into any liquid) and sharp edges. Do not use the appliance with wet hands.
- Do not operate the machine without supervision.
- Do not use the appliance if it is damaged.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- The device is used for massaging the skin of adults.
- The device should only be used for this purpose.
- Do not use the device at eye level or on other sensitive parts of the body.
- The massager should not be used under a blanket or in other poorly ventilated places.

Installing the Battery

The battery compartment is located at the bottom end of the handle.

- Release the screw on the cover.
- Insert 2 batteries of type AA • UM-3 • R6.
- Please note the polarity! See the information on the device.

Batteries (not supplied)

- If the device is not used for a long period of time, remove the batteries to avoid the spilling of battery acid.
- Different battery types or new and used batteries must not be used together.

WARNING:

Do not expose the batteries to high temperatures or direct sunlight. Never throw batteries into the fire. There is a danger of explosion!

Function switch

- 0 = off
- 1 = gentle massage
- 2 = more intensive massage

Cleaning

CAUTION:

- Do not use abrasive detergents.
- The machine must never be immersed in water!
- Clean the machine only with a slightly damp cloth.

Net weight:.....0,19 kg

Battery Disposal

Your massage device includes batteries. When these batteries have worn out, please comply with the following:

CAUTION:

- Standard/rechargeable batteries must not be disposed of in the normal domestic waste!
- As a consumer you are legally obliged to return used standard/rechargeable batteries.
- The batteries have to be removed from the device.
- You can return your standard/rechargeable batteries free of charge wherever the batteries were bought.
- Information on public collection points can be obtained from your local authority.

Environmental Protection

When you wish to dispose of the device, please do not put it in the standard domestic waste. Instead you should take it to an official collection point for recycling. In this way you help to protect the environment.

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